



**Restriction Map and Multiple Cloning Site (MCS) of pHIS2 Vector.** Unique restriction sites are in bold.

**\*Note:** The *Sac II* site is resistant to cleavage.

### Description

pHIS2 is a reporter vector that can be used in yeast one-hybrid assays to identify and characterize DNA-binding proteins. The vector was specifically designed for use with the BD Matchmaker™ One-Hybrid Library Construction & Screening Kit (#K1617-1). It contains a *HIS3* nutritional reporter gene, located downstream of a multiple cloning site (MCS) and the minimal promoter of the *HIS3* locus ( $P_{\min HIS3}$ ). *Cis*-acting DNA sequences, or DNA target elements, can be inserted into the MCS and used as baits to screen GAL4 AD/cDNA fusion libraries for proteins that interact with the target sequence. A protein-DNA (or one-hybrid) interaction can be detected by performing the assay in a yeast strain such as Y187 that is auxotrophic for histidine. Positive one-hybrid interactions drive expression of the *HIS3* reporter gene, which enables the host cell to grow on histidine-deficient media.

营养缺陷型

In the absence of activation, the constitutive *HIS3* expression from  $P_{\min HIS3}$  is very low. During library screening, the leaky expression of *HIS3* is controlled by adding 3-amino-1,2,4-triazole (3-AT) to the medium. The concentration of 3-AT needed to fully suppress leaky *HIS3* expression must be determined empirically for each DNA target element.

pHIS2 can be maintained in both yeast and bacteria. It contains an autonomous replication sequence (*ARS4*) and *TRP1* nutritional marker for replication and selection in yeast (1, 2); it contains a Col E1 origin and a kanamycin resistance gene ( $Kan^r$ ) for propagation and selection in *E. coli*. The centromeric sequence *CEN6* ensures proper segregation of the plasmid during cell division in yeast (1, 2).

## Use

To use pHIS2 in a one-hybrid assay, clone one or more copies of a *cis*-acting DNA target sequence into the MCS. Then introduce the plasmid into competent yeast cells using the transformation protocol in the BD Matchmaker Library Construction & Screening Kits User Manual (PT3529-1). In contrast to the original BD Matchmaker One-Hybrid System, this reporter vector does not need to be integrated into the yeast genome. Instead, it is maintained as an episome throughout the assay. Inserting your target element may alter the level of background *HIS3* expression. Therefore, constructs should be tested for background (leaky) *HIS3* expression before you start a one-hybrid analysis. Background growth due to leaky *HIS3* expression is controlled by adding 3-AT to the selection medium, as described in the User Manual (PT3529-1).

## Location of Features

- Multiple cloning sites: 1–41
- *HIS3* gene: 152–811
- Fragment containing the *HIS3* 3' UTR & Termination sequence: 812–1446
- *TRP1* gene: 2855–3529
- Fragment containing the Col E1 *E. coli* origin of replication: 4165–4615
- Kanamycin-resistance gene: 5605–4811
- *CEN6/ARS4* sequences: 6254–5737

## Propagation in *E. coli*

- Suitable host strains: DH5 $\alpha$  and other general purpose strains.
- Selectable marker: plasmid confers resistance to kanamycin (50  $\mu$ g/ml) to *E. coli* hosts.
- *E. coli* replication origin: Col E1
- Copy number: low

## Propagation in *S. cerevisiae*

- Suitable host strains: AH109(MAT $\alpha$ ), Y187(MAT $\alpha$ ), Y190(MAT $\alpha$ ), SFY526(MAT $\alpha$ ), CG1945(MAT $\alpha$ ), HF7c(MAT $\alpha$ )
- Selectable marker: *TRP1*
- *S. cerevisiae* origin: *CEN6/ARS4*

## References

1. Sikorski, R. S. & Hieter, P. (1989) *Genetics* **122**:19–27.
2. Rose, M. D. & Broach, J. R. (1991) *Methods Enzymol.* **194**:195–230.

**Note:** The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by BD Biosciences Clontech. This vector has not been completely sequenced.

## Notice to Purchaser

This product is intended for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. BD Biosciences Clontech products may not be resold, modified for resale, or used to manufacture commercial products without written approval of BD Biosciences Clontech.

BD and BD Logo are trademarks of Becton, Dickinson and Company. ©2003 BD